

أحمر LED عصا العلاج المضادة للشيخوخة، ضوء **LED** نانومتر 660 نانومتر ضوء 850

سبائك الألومنيوم مع اللون الأسود المطلي

القاعدة : ه 27

التبريد التلقائي مع درع القالب

أداء ثابت لرقائق الأشعة تحت الحمراء

.شريحة 3 واط لكل جهاز كمبيوتر

قطعة من الأشعة تحت الحمراء: 850 نانومتر 9

قطعة من الضوء الأحمر: 660 نانومتر 9














54W

Powerful
Red LED Light



Infrared light enhance cell vitality

Under infrared irradiation, capillary dilation, metabolism and cell viability and regeneration ability were improved




54 watt infrared light source

Because the distribution of heat on the irradiation area is more uniform, it covers a wider area, which can effectively relieve pain



Perfect Wavelength Combination for Light Therapy

PCS high-power deep red 660nm and 9PCS high-power deep red 850nm(looks dark, but more efficient) can penetrate deeper into muscles and body tissues, so as to more effectively relieve pain



Stimulate the circulation of blood

The infrared light penetrates into the deep layer of the skin, promotes the thermal cycle, increases the phagocytosis of white blood cells, and improves the immunity of the body